

## Prescreen/Rescreen workflow on an ACQUIFER Imaging Machine



Experiment Default Channels 2 X Positions 24W Dimension 2D Autofocus None Run Completed Smart Imaging Disable System Settings... Plate Viewer Launch App



## Objective

## Camera Binning

2 X 4 X 10 X 20 X

1 x 1


## Channels

	On / Off	Channel	Exp	Exposure [ms]	Illumination Power [%]	Modify
1	<input checked="" type="checkbox"/> On	BF	<input checked="" type="checkbox"/>	20	30	Edit Remove


+ Add

## Well Navigation

Go



## Objective Navigation



Y

Z

X

Steps

XY 3.000 mm

Z 250.0 μm

Z-Stack Center

Go To Z-Stack Center

Current Position

X 113.157 mm

Y 20.326 mm


Z 21500.1 μm

## Histogram

Manual Brightness

1314

16843



## Autofocus

No data to plot

# Sequential Prescreen/Rescreen workflow - python



```

%% Import and open tcpip communication
from acquifer import tcpip, scripts
from ScriptUtils import PixelPosition
import MTM, cv2, os

%% Checks
hasUpdatedValues = input("Did you update values for path_prescreen, path_rescreen,
if not hasUpdatedValues.strip().lower() in ["y", "yes"]:
    raise InterruptedError("Update values first")

%% Input scripts
path_prescreen = r"C:\Users\Administrator\Downloads\2X-script.imsf"
path_rescreen = r"C:\Users\Administrator\Downloads\4X-script.imsf"
path_template = r"C:\Users\Administrator\Downloads\medaka_crop.tif"

zref = 21500.1

%% Run prescript
scope = tcpip.TcpIp()
directory_prescreen = scope.runScript(path_prescreen)
    
```

```

listPositions = []
for filename in os.listdir(directory_prescreen):
    if not filename.endswith(".tif"):
        continue

    if not "C01" in filename:
        continue

    filepath = os.path.join(directory_prescreen, filename)

    image = cv2.imread(filepath, -1)

    hits = MTM.matchTemplates(listTemplates=[("template", template)],
                                image = image,
                                N_object = 1,
                                score_threshold=0.5,
                                method=cv2.TM_CCOEFF_NORMED,
                                maxOverlap=0)

    if len(hits) == 0:
        continue

    print(hits)

    score = hits["Score"][0]
    if score < 0.5:
        continue

    bbox = hits["BBBox"][0]
    x,y,width, height = bbox
    bboxCenter_x = int(x + width/2)
    bboxCenter_y = int(y + height/2)

    # Crop detected region and save it
    foundImage = image[x : x+width, y : y+height]
    cv2.imwrite(os.path.join(directory_detected, filename), foundImage)

    # Create a pixel position and add it to the list
    position = PixelPosition(bboxCenter_x, bboxCenter_y, float(zref), filepath)
    listPositions.append(position)

%% Update positions and run script with new positions
script = scripts.replacePositionsInScriptFile(path_rescreen, listPositions)

scope = tcpip.TcpIp()
scope.runScript(script)
    
```

www.BANDICAM.COM

ACQUIFER

Experiment Default Channels 2 X Positions 24W Dimension 2D Autofocus None Run Completed Smart Imaging Disable System Settings... Plate Viewer Launch App

Objective: 2 X, 4 X, 10 X, 20 X. Camera Binning: 1 x 1

Channels:

	On / Off	Channel	Exp	Exposure [ms]	Illumination Power [%]	Modify
1	On	BF	✓	20	30	Edit Remove

Add

Well Navigation

Navigation

Steps: XY 3.000 mm, Z 250.0 µm

Current Position: X 113.157 mm, Y 20.326 mm, Z 21500.1 µm

Manual Brightness: 16799

Autofocus

Bandicam Settings:

Screen Recording Mode - Record Desktop screen and Web video.

Record:

- ☒ Record/Stop Hotkey: F12
- ☐ Pause Hotkey: Shift+F12
- ☒ Show mouse cursor
- ☒ Add mouse click effects
- ☐ Add webcam overlay

Format - AVI:

Video: H264 - NVIDIA® NVENC (VBR) (Auto), Full Size, 60.00fps, 80q

Audio: MP3 - MPEG-1 Layer 3, 48.0KHz, stereo, 192kbps

Remove the watermark and 10 minutes recording limit

LUT: None, Mouse-Position: X: 1235 Y: 1333, Grey-Value: 13115, Zoom: 0.56, Width: 2048 px, Height: 2048 px

3:53 PM 1/18/2024

## Pros

- simple / high-level
- single short python file (+ protocol for control software)

## Cons

- Sequential : no „real time“ feedback
- No direct access to the microscope parameters (exposure, light-source intensity...)

Note : „low-level“/fine grain approach possible too but more complex/verbose

# ACQUIFER Imaging Machine

Various options for smart-microscopy

